

Bondller

Double Sided Tissue Tape 2501

Product Description

Superior adhesive strength, permanent mount, can be curved, repulsion. High adhesion to wide range of materials(like plastic, metal...) Good season resistant, low and high sevice temperautre ability

Technical Data

| Test standard | China GB | ASTM | FINAT | JIS | |
|--------------------|-------------|-----------|-----------|-----------|--|
| Tape Thickness(mm) | 0.14 | | | | |
| Color | Translucent | | | | |
| Adhesion Power | 15N/25mm | 15N/25mm | 15N/25mm | 15N/25mm | |
| Initial Tack | 13#(ball) | 3cm | 15N/25mm | 13#(ball) | |
| Holding Power | >72 hours | >72 hours | >72 hours | >72 hours | |

Application

Bonding punched metal or plastic nameplates, fixing plastic display panel(like automobile meter panels and air-conditioners), and also permanent laminating household electric appliances components and parts.



Construction



Modified solvent base acrylic adhesive Tissue 14g/m2 Modified solvent base acrylic adhesive PE paper release liner 115g/m2(with print LUXKING)

Temperature

Application (Minimum): 10° Service (long term): -20° to 60° Short term: 80°

Standard Size

| | Tape width | Usable width | Length | | |
|-------------|--------------|--------------|-----------------|--|--|
| Jumbo rolls | 1040, 1260mm | 1020, 1240mm | 500, 700, 1000M | | |
| Log rolls | 1040, 1260mm | 1020, 1240mm | 50, 100M | | |

Surface Preparation

It is essential, as with all pressure sensitive tapes, that the surface to which the tape is applied be clean, dry and free of grease and oil

Storage and Shelf Life The shelf life of the products is one years when stored at 18-25°C, 30-70% releative humidity and out of direct sunlight, in good packaging.